

CZR DENTAL PORCELAIN for Zirconia Framework Restorations

CZR BASIC FULL KIT / CZR-SS BASIC FULL KIT

CE 0120

CZR BASIC FULL KIT is an excellent kit composed of basic colors for all ceramic crowns and bridges, using Zirconia frameworks.

CZR-SS BASIC FULL KIT has the same shade selection as CZR BASIC FULL KIT. Shade Base Stains are contained instead of Shade Base porcelains.

Coefficient of thermal expansion
(50–500°C 10⁻⁶ K⁻¹)

EX-3	Cerabien	CZR
12.4	6.8	9.1

SB (powder) / SS (stain)	10g each (SS 3g each)	18 colors
BODY		18 colors
ENAMEL		4 colors
TRANSLUCENT		4 colors
CERVICAL		4 colors
MODIFIER		13 colors
ADD-ON		2 colors
TECHNICAL INSTRUCTIONS		1 copy



CZR BASIC FULL KIT

CZR MARGIN KIT

CE 0120

In addition to ordinary usage, MARGIN PORCELAIN can be used to correct a chip in the Zirconia framework. The use of MAGIC FORMER as forming liquid is highly recommended.

BASIC SHADE	10g each	14 colors
SPECIAL COLOR		3 colors
DILUTION POWDER (MDL)		1 color
RETOUCHING POWDER (MRP)		1 color



CZR OPACIOUS BODY COMPLETE KIT

CE 0120

OPACIOUS BODY can be used to achieve delicate colors and to control translucency.

BASIC SHADE	10g each	16 colors
SPECIAL COLOR		OB White OB Orange OB Pale Pink OBE



CZR LUSTER CCV KIT

CE 0120

LUSTER PORCELAIN reproduces natural enamel in both surface features and color.

- LUSTER PORCELAIN reproduces the fine surface structure and luster of natural teeth.
- A unique combination of fine surface particles produces a selective reflection of light that results in the same opalescence seen in natural teeth.

LT ₀	LT ₁	10g each
T Blue	AQUA BLUE 1	
AQUA BLUE 2	LT NATURAL	
LT SUPER GRAY	CREAMY ENAMEL	
SUN BRIGHT	INCISAL AUREOLA	
CREAMY WHITE	LY YELLOW	
CCV-1	CCV-2	
CCV-3	CCV-4	



CZR DENTAL PORCELAIN for Zirconia Framework Restorations

CZR ESTHETIC WHITE KIT / CZR-SS ESTHETIC WHITE KIT CE 0120

CZR ESTHETIC WHITE KIT is for creating whiter shades rather than NW shades. It is not only pure bleaching shades but also natural and beautiful white shades.

CZR-SS ESTHETIC WHITE KIT has a Shade Base Stain (SSB₁) instead of Shade Base SBB₁.

SB (powder) / SS (stain)	10g each (SS 3g each)	SB White, SBB ₁ / SS White, SSB ₁
BODY		EW ₀₀ B, EW ₀ B, EWB, EWYB
ENAMEL		Silky E ₁ , Silky E ₂
LUSTER TRANSLUCENT		ELT ₁ , ELT ₂ , ELT ₃



CZR-SS ESTHETIC WHITE KIT

CZR TISSUE KIT CE 0120

CZR TISSUE KIT is complete kit of 7 Noritake Tissue colors. All colors are very useful for implant work.

TISSUE	10g each	Tissue1~7
--------	----------	-----------



CZR INTERNAL STAIN COMPLETE KIT CE 0120

Using INTERNAL STAIN, it is possible to replicate many of the most intricate three-dimensional patterns of shade variation found in natural teeth. IS FLUORO is added to the kit.

WHITE	CERVICAL 2	GRAY	3g each
INCISAL BLUE 1	CERVICAL 3	FLUORO	
INCISAL BLUE 2	REDDISH BROWN	A+	
MAMELON ORANGE 1	REDDISH BROWN	B+	
MAMELON ORANGE 2	SALMON PINK	C+	
CERVICAL 1	RED	D+	
BRIGHT (DILUTION)		3g	1 color
IS LIQUID		10ml	1 bottle
IS COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy



CZR EXTERNAL STAIN KIT CE 0120

EXTERNAL STAIN reproduces the great variety of color shades observed on the surface of natural teeth.

PURE WHITE	ORANGE 1	SALMON PINK	3g each
GRAY	ORANGE 2	PINK	
BLACK	CERVICAL 1	RED	
BLUE	CERVICAL 2	A+	
GREEN 1	CERVICAL 3	B+	
GREEN 2	EARTH BROWN	C+	
YELLOW	REDDISH BROWN	D+	
GLAZE		10g	1 color
ES LIQUID		10ml	1 bottle
ES COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy

